

## Atlantis® (IO) FLO-S compatibility chart

The scan body Atlantis (IO) FLO-S is used when scanning for Atlantis suprastructures. It can be used on implant/abutment level for all positions in the mouth of the patient or on the master model during scanning. Refer to instructions for use for further information, 2663-IFU Atlantis IO FLO-S.

This chart describes what implants/abutments the Atlantis (IO) FLO-S is compatible with and the content of Dentsply Sirona Implants Atlantis suprastructures (IO) FLO-S libraries. The libraries are used for digital design guidance for Atlantis suprastructures and for building and printing models. The libraries can be used for 3Shape and Exocad design software.

**Note:** The FLO-S scan body without the IO marking can only be used for model scanning and is not approved for intraoral scanning. The IO FLO-S scan body will be used for both intraoral and model scanning.

Each Atlantis (IO) FLO-S reference number is sold as a kit containing a tray with 6 (IO) FLO-S for storing the non-sterile (IO) FLO-S and a Tag (to be placed in the Atlantis FLO-S case (REF 36701) for 1 specific connection.



Manufacturer	Implant system	FLO-S id	FLO-S Movex no (Kit)	Intraoral scanning	Library information			
					Available in library	Elos Accurate Analog for Printed Models	Elos Accurate Hybrid Base Non-Engaging	Elos Accurate Hybrid Base H Non-Engaging
Dentsply Sirona	<b>IMPLANT LEVEL</b>							
	<b>EV Connection: PrimeTaper EV Implant / OmniTaper EV Implant / Astra Tech Implant EV</b>							
	EV Connection (XS) Green	A04K	68020018	YES	YES	PMA-ATE30-1 (68020083)		
	EV Connection (S) Purple	A04L	68020019	YES	YES	PMA-ATE36-1 (68020084)		
	EV Connection (M) Yellow	A04M	68020020	YES	YES	PMA-ATE42-1 (68020085)		
	EV Connection (L) Blue	A04N	68020021	YES	YES	PMA-ATE48-1 (68020086)		
	EV Connection (XL) Brown*	A04O	68020022	YES	YES	PMA-ATE54-1 (68020087)		
	<b>Astra Tech Implant EV Profile</b>							
	4.2 Yellow	A04M	68020020	YES	YES	PMA-ATE42P-1 (68020088)		
	4.8 Blue	A04N	68020021	YES	YES	PMA-ATE48P-1 (68020089)		
	<b>Astra Tech Implant System TX</b>							
	OsseoSpeed TX 3.0 – Yellow	A04C	68020023	YES	YES	PMA-ATO30-1 (68020090)		
	OsseoSpeed TX 3.5, 4.0 – Aqua	A04D	68020024	YES	YES	PMA-ATO35-1 (68020091)		
	OsseoSpeed TX 4.5, 5.0 – Lilac	A04E	68020025	YES	YES	PMA-ATO45-1 (68020092)		
	OsseoSpeed TX Profile 4.5/5.0	A04E	36741-6		YES			
	<b>Xive / Frialit</b>							
	Xive S 3.0	F01A	68020026	YES	YES	PMA-DXI30-1 (68020093)		
	Xive S 3.4 / Frialit 3.4	F01B	68020027	YES	YES	PMA-DXI34-1 (68020094)		
	Xive S 3.8 / Frialit 3.8	F01C	68020028	YES	YES	PMA-DXI38-1 (68020095)		
	Xive S 4.5 / Frialit 4.5	F01D	68020029	YES	YES	PMA-DXI45-1 (68020096)		
	Xive S 5.5 / Frialit 5.5	F01E	68020030	YES	YES	PMA-DXI55-1 (68020097)		

\* Note: Not applicable for PrimeTaper EV Implant.

All trademarks and company names are the property of their respective owners.

Elos Accurate® is a registered trademark.

**Note:** Atlantis (IO) FLO-S assortment can vary from country to country. Please check availability with your local Customer Service Representative.

Please turn over →

Manufacturer	Implant system	FLO-S id	FLO-S Movex no (Kit)	Intraoral scanning	Library information			
					Available in library	Elos Accurate Analog for Printed Models	Elos Accurate Hybrid Base Non-Engaging	Elos Accurate Hybrid Base H Non-Engaging
Dentsply Sirona	<b>ABUTMENT LEVEL</b>							
	<b>EV Connection: DS PrimeTaper Implant / DS OmniTaper Implant / Astra Tech Implant EV</b>							
	MultiBase Abutment EV	A04S	68020034	YES	YES	PMA-MUA45-1 (68020099)	HBN-MUA45-1	HBN-MUA45H-1
	Uni Abutment EV	A04F	68020031	YES	YES	PMA-AUE33-1 (68020100)	HBN-AUE33-1	HBN-AUE33H-1
	<b>Astra Tech Implant System TX</b>							
	20° UniAbutment	A04A	68020032	YES	YES	PMA-ATU20-1 (68020101)	HBN-ATU20-1	HBN-ATU20H-1
	45° UniAbutment	A04B	68020033	YES	YES	PMA-ATU45-1 (68020102)	HBN-ATU45-1	
	<b>Xive</b>							
	Xive MP/TG 3.4,3.8	F01F	36745-6		YES			
	Xive MP/TG 4.5	F01G	36746-6		YES			
	Xive MP 5.5	F01H	36747-6		YES			
	<b>Ankylos</b>							
	Ankylos Balance Base Narrow D4.2	A03A	68020035	YES	YES	PMA-DBA40-1 (68020103)	HBN-DBA40-1	HBN-DBA40H-1
	Ankylos Balance Base D5.5	A03B	36730-6		YES			
MIS	<b>IMPLANT LEVEL</b>							
	Internal Hex Narrow SEVEN – M4	M04D	68020043	YES	YES			
	Internal Hex Standard SEVEN – M4	M04A	68020044	YES	YES			
	Internal Hex Wide SEVEN – M4	M04B	68020045	YES	YES			
	Conical Connection Narrow C1	M04E	68020046	YES	YES			
	Conical Connection Standard C1 – V3	M04F	68020047	YES	YES			
	Conical Connection Narrow V3	M04E	68020046	YES	YES			
	Conical Connection Wide C1	M04G	68020048	YES	YES			
	<b>ABUTMENT LEVEL</b>							
	Multi-Unit Abutment	M04H	68020049	YES	YES			
	CONNECT abutment D4	M04I	68021007					
	CONNECT abutment D5.7	M04J	68021008					
	BioHorizons	<b>IMPLANT LEVEL</b>						
Internal/ Tapered 3.5		B03A	36710-6					
Internal/ Tapered 4.5		B03B	36711-6					
Internal/ Tapered 5.7		B03C	36712-6					
<b>ABUTMENT LEVEL</b>								
Multi-unit abutment		A04S	68020034					

Please turn over 

All trademarks and company names are the property of their respective owners.

Elos Accurate® is a registered trademark.

**Note:** Atlantis (IO) FLO-S assortment can vary from country to country. Please check availability with your local Customer Service Representative.

Manufacturer	Implant system	FLO-S id	FLO-S Movex no (Kit)	Intraoral scanning	Library information			
					Available in library	Elos Accurate Analog for Printed Models	Elos Accurate Hybrid Base Non-Engaging	Elos Accurate Hybrid Base H Non-Engaging
Biomet 3i	<b>IMPLANT LEVEL</b>							
	Certain 3.25, 4/3	301F	36716-6					
	Certain 4.0, 5/4	301H	36717-6					
	Certain 5.0	301J	36718-6					
	Certain 6.0	301L	36719-6					
	Certain Prevail 3/4/3, 4/3	301F	36716-6					
	Certain Prevail 4/5/4, 5/4	301H	36717-6					
	Certain Prevail 5/6/5, 6/5	301J	36718-6					
	Certain XP4/5	301J	36718-6					
	Certain XP 5/6	301L	36719-6					
	MicroMiniplant 3.25 / Miniplant 3.25	301G	36721-6					
	Osseotite 3.75 / 4.0	N01F	68020009					
	Osseotite 5.0, XP 4/5	301K	36723-6					
	Osseotite 6.0, XP 5/6	301M	36724-6					
<b>ABUTMENT LEVEL</b>								
Low Profile Abutment	301E	36720-6						
Camlog	<b>IMPLANT LEVEL</b>							
	Screw-line Implant 3.3	C01E	36725-6					
	Screw-line/Root-line Implant 3.8	C01G	36726-6					
	Screw-line/Root-line Implant 4.3	C01H	36727-6					
	Screw-line/Root-line Implant 5.0	C01I	36728-6					
	Screw-line/Root-line Implant 6.0	C01K	36729-6					
Keystone Dental	<b>IMPLANT LEVEL</b>							
	PrimaConnex SD 3.3,3.5	K01A	36755-6					
	PrimaConnex RD 4.0,4.1	K01B	36754-6					
	PrimaConnex WD 5.0	K01C	36756-6					
	Genesis 3.8	K01A	36755-6					
	Genesis 4.5	K01B	36754-6					
	Genesis 5.5, 6.5	K01C	36756-6					
Megagen	<b>IMPLANT LEVEL</b>							
	AnyOne 3.5	M03A	68021002					
	AnyOne 4.0, 4.5, 5.0, 5.5, 6.0, 6.5, 7.0, 7.5, 8.0	M03A	68021002					
	<b>ABUTMENT LEVEL</b>							
	Multi-unit Abutment (N_Type)	A04S	68020034					
	Multi-unit Abutment (S_Type)	M03G	68021004					
	Octa Abutment 3.8	M03H	68021005					
	Octa Abutment 4.8	S05C	36772-6					
	Octa Abutment 5.8	M03I	68021006					

Please turn over 

All trademarks and company names are the property of their respective owners.

Elos Accurate® is a registered trademark.

**Note:** Atlantis (IO) FLO-S assortment can vary from country to country.  
Please check availability with your local Customer Service Representative.

Manufacturer	Implant system	FLO-S id	FLO-S Movex no (Kit)	Intraoral scanning	Library information			
					Available in library	Elos Accurate Analog for Printed Models	Elos Accurate Hybrid Base Non-Engaging	Elos Accurate Hybrid Base H Non-Engaging
<b>Neodent</b>	<b>IMPLANT LEVEL</b>							
	GM Implant 3.5, 3.75, 4.0, 4.3	N03A	68021009					
	GM Implant 5.0, 6.0, 7.0	N03A	68021009					
	<b>ABUTMENT LEVEL</b>							
	GM Abutment	N03C	68021010					
	GM Micro Abutment	N03B	68021011					
GM Mini Conical Abutment	N01H	36759-6						
<b>Nobel Biocare</b>	<b>IMPLANT LEVEL</b>							
	NobelActive NP 3.5	N01A	68020050	YES	YES	PMA-NBA35-1 (68020119)		
	NobelActive RP 4.3, 5.0	N01B	68020051	YES	YES	PMA-NBA43-1 (68020120)		
	NobelActive WP 5.5	N01O	68020052	YES	YES	PMA-NBA60-1 (68020121)		
	NobelReplace NP 3.5	N01J	68020053	YES	YES	PMA-NBR35-1 (68020122)		
	NobelReplace RP 4.3	N01K	68020054	YES	YES	PMA-NBR43-1 (68020123)		
	NobelReplace WP 5.0	N01L	68020055	YES	YES	PMA-NBR50-1 (68020124)		
	NobelReplace 6.0	N01M	68020056	YES	YES	PMA-NBR60-1 (68020125)		
	Branemark System NP 3.3	N01E	68020057	YES	YES	PMA-BRA35-1 (68020126)		
	Branemark System RP 3.75, 4.0	N01F	68020009	YES	YES	PMA-BRA41-1 (68020127)		
	<b>ABUTMENT LEVEL</b>							
	Multi-Unit Abutment RP	N01H	68020058	YES	YES	PMA-MUA45-1 (68020099)		
	<b>Straumann</b>	<b>IMPLANT LEVEL</b>						
Bone Level 3.3 NC		S05G	68020059	YES	YES	PMA-SBO33-1 (68020129)		
Bone Level 4.1, 4.8 RC		S05H	68020060	YES	YES	PMA-SBO41-1 (68020130)		
Standard/Standard Plus 4.8 RN (Regular Neck)		S05D	68020061	YES	YES	PMA-SSY48-1 (68020131)		
Standard/Standard Plus 4.8 WN (Wide Neck)		S05F	68020062	YES	YES	PMA-SSY65-1 (68020132)		
BLX Implant RB 3.75, 4.5		S05S	68021015					
BLX Implant WB 5.5, 6.5		S05S	68021015					
<b>ABUTMENT LEVEL</b>								
Monotype 8° Conus		S05A	36771-6		YES			
RN Abutment Level		S05C	36772-6		YES			
WN Abutment Level		S05E	36777-6		YES			
Bone Level Multi-Base Abutment D3.5		S05I	36768-6		YES			
Bone Level Multi-Base Abutment D4.5		S05J	36769-6		YES			
Bone Level Multi-Base Angled Abutment		S05K	36770-6		YES			
Screw-Retained Abutment 3.5		S05O	68020069	YES	YES	PMA-SSA35-1 (68020135)		
Screw-Retained Abutment 4.6	S05P	68020070	YES	YES	PMA-SSA46-1 (68020136)			

Please turn over 

All trademarks and company names are the property of their respective owners.

Elos Accurate® is a registered trademark.

**Note:** Atlantis (IO) FLO-S assortment can vary from country to country. Please check availability with your local Customer Service Representative.

Manufacturer	Implant system	FLO-S id	FLO-S Movex no (Kit)	Intraoral scanning	Library information			
					Available in library	Elos Accurate Analog for Printed Models	Elos Accurate Hybrid Base Non-Engaging	Elos Accurate Hybrid Base H Non-Engaging
Zimmer Dental	<b>IMPLANT LEVEL</b>							
	Tapered Screw-Vent™ 3.5 / Screw Vent™ 3.3, 3.7	B03A	36710-6					
	Tapered Screw-Vent™ 4.5 / Screw Vent™ 4.5	B03B	36711-6					
	Tapered Screw-Vent™ 5.7	B03C	36712-6					
	<b>ABUTMENT LEVEL</b>							
	Tapered Abutment	Z01E	36779-6					
Spectra Cone Abutment	Z01F	36778-6						

All trademarks and company names are the property of their respective owners.

Elos Accurate® is a registered trademark.

**Note:** Atlantis (IO) FLO-S assortment can vary from country to country. Please check availability with your local Customer Service Representative.

Tool for Elos Accurate® Analog for Printed Models	REF
Elos Accurate Analog Pliers	PMA-AIT-1 (68020141)
Elos Accurate Analog Insertion Screw	PMA-AIS-5 (68020142)
Elos Accurate Analog Insertion Pin	PMA-AIP-2 (68020143)

